

Rabbit Anti-Mouse Cd45 Monoclonal Antibody For If/Icc And Fc - P06800

Item Number: CM0008735

Introduction

-conjugated rabbit monoclonal antibody against mouse CD45 (PTPRC), validated for immunofluorescence (IF/ICC) and flow cytometry (FC). Detects leukocyte common antigen (L-CA) in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD45 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GenelD (Human)	5788
GenelD	19264
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 26-566 of mouse CD45 (NP_001104786.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag