

Rabbit Anti-Mouse Cd45.1 Mab - P06800

Item Number: CM0005600

Introduction

-conjugated rabbit monoclonal antibody against mouse CD45.1 (P06800). Validated for flow cytometry, this antibody enables sensitive detection of leukocyte populations in mouse models. Ideal for immunophenotyping and T-cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD45.1 mAb
Remarks/Alias	loc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 311-410 of mouse CD45.1 (NP_001104786.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag