

Rabbit Anti-Mouse Cd45.1 Monoclonal Antibody For Fc - P06800

Item Number: CM0007728

Introduction

Rabbit monoclonal antibody against mouse CD45.1, validated for flow cytometry (FC). Ideal for immunophenotyping, T-cell activation studies, and distinguishing CD45 allotypes in mouse models.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD45.1 mAb
Remarks/Alias	Ioc; B220; Cd45; L-CA; Ly-5; T200; CD45R; Lyt-4
Species	Mouse
GeneID (Human)	5788
GeneID	19264
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130 kDa/140 kDa/144 kDa
Shipping	Ice bag