

Rabbit Anti-Mouse Cd45.1 Monoclonal Antibody For Flow Cytometry - P06800

Item Number: CM0011834

Introduction

High-quality -conjugated rabbit anti-mouse CD45.1 monoclonal antibody for flow cytometry, targeting the leukocyte common antigen. Validated in mouse samples, suitable for immunophenotyping and T-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD45.1 mAb
Remarks/Alias	loc, B220, Cd45, L-CA, Ly-5, T200, CD45R, Lyt-4
Species	Mouse
Gene ID (Human)	5788
Gene ID	19264
Immunogen	Synthetic peptide, A synthetic peptide corresponding to a sequence within amino acids 311-410 of mouse CD45.1 (NP_001104786.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P06800
Protein Weight	130kDa/140kDa/144kDa
Shipping	Ice bag