

Rabbit Anti-Mouse Cd47 Mab For Flow Cytometry - Q61735

Item Number: CM0009440

Introduction

High-specificity -conjugated rabbit monoclonal antibody against mouse CD47 (IAP).
Recombinant immunogen (aa 19-135). Reacts specifically with mouse CD47. Designed for flow cytometry (FC) in immunology and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD47 mAb
Remarks/Alias	IAP; Itgp; 9130415E20Rik; B430305P08Rik
Species	Mouse
GeneID (human)	961
GeneID	16423
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of mouse CD47 (NP_034711.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61735
Protein Weight	33KDa/35KDa
Shipping	Ice bag