

Rabbit Anti-Mouse Cd48 Monoclonal Antibody For Fc - P18181

Item Number: CM0008000



Introduction

-conjugated rabbit monoclonal antibody against mouse CD48 (SLAMF2/BCM1) for flow cytometry. Validated for mouse reactivity; target protein weight 27kDa. Suitable for T-cell and NK-cell immunology studies.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Mouse CD48 mAb
Remarks/Alias	BCM1; BLAST; Bcm-1; BLAST1; SLAMF2; Sgp-60; BLAST-1; MEM-102
Species	Mouse
GeneID (Human)	962
GeneID	12506
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 23-216 of mouse CD48 (NP_031675.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18181
Protein Weight	27kDa
Shipping	Ice bag