

Rabbit Anti-Mouse Cd48 Monoclonal Antibody For Fc - P18181

Item Number: CM0007875



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD48 for flow cytometry (FC). Immunogen: recombinant mouse CD48 extracellular domain (aa 23-216). CD48, a GPI-anchored SLAM family receptor, mediates immune cell activation and NK cell cytotoxicity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD48 mAb
Remarks/Alias	BCM1; BLAST; Bcm-1; BLAST1; SLAMF2; Sgp-60; BLAST-1; MEM-102
Species	Mouse
Gene ID (Human)	962
Gene ID	12506
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-216 of mouse CD48 (NP_031675.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18181
Protein Weight	27kDa
Shipping	Ice bag