

Rabbit Anti-Mouse Cd48 Monoclonal Antibody For Flow Cytometry - P18181

Item Number: CM0005034



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD48 (SLAMF2). Validated for flow cytometry (FC). Ideal for immunology research involving T-cell and NK-cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD48 mAb
Remarks/Alias	BCM1; BLAST; Bcm-1; BLAST1; SLAMF2; Sgp-60; BLAST-1; MEM-102
Species	Mouse
GeneID (human)	962
GeneID	12506
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 23-216 of mouse CD48 (NP_031675.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18181
Protein Weight	27kDa
Shipping	Ice bag