

Rabbit Anti-Mouse Cd5 Monoclonal Antibody For Flow Cytometry - P13379

Item Number: CM0009190

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD5 (Ly-1). Validated for flow cytometry, detects a 54 kDa T-cell and B1a cell marker. Suitable for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD5 mAb
Remarks/Alias	Ly-1; Ly-A; Ly-12; Lyt-1
Species	Mouse
GeneID (Human)	921
GeneID	12507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of mouse CD5 (NP_031676.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13379
Protein Weight	54kDa
Shipping	Ice bag