

Rabbit Anti-Mouse Cd5 Monoclonal Antibody For Flow Cytometry - P13379

Item Number: CM0009104

Introduction

Rabbit anti-Mouse CD5 Monoclonal Antibody validated for Flow Cytometry (FC) on Mouse samples. Targets T-cell surface glycoprotein CD5 (Ly-1). Ideal for immunology research on T/B cell signaling and homeostasis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD5 mAb
Remarks/Alias	Ly-1; Ly-A; Ly-12; Lyt-1
Species	Mouse
GeneID (Human)	921
GeneID	12507
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of mouse CD5 (NP_031676.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13379
Protein Weight	54kDa
Shipping	Ice bag