

Rabbit Anti-Mouse Cd56/Ncam Monoclonal Antibody For Fc - P13595

Item Number: CM0009630



Introduction

This -conjugated rabbit monoclonal antibody specifically recognizes mouse CD56/NCAM, making it an ideal tool for flow cytometric analysis of neural cell adhesion molecule expression. It detects multiple isoforms with molecular weights of 67, 80, 93, and 119 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD56/NCAM mAb
Remarks/Alias	CD56; Ncam; E-NCAM; NCAM-1
Species	Mouse
GeneID (human)	4684
GeneID	17967
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-711 of mouse CD56 (NP_001106675.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P13595
Protein Weight	67kDa/80kDa/93kDa/119kDa
Shipping	Ice bag