

Rabbit Anti-Mouse Cd59A Monoclonal Antibody For Fc - 055186

Item Number: CM0011951

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD59a for flow cytometry. Detects a 14kDa GPI-anchored MAC inhibitor; cross-reacts with mouse. Ideal for complement research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD59a mAb
Remarks/Alias	Cd59; MACIF; MAC-IP
Species	Mouse
GeneID	12509 (mouse)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-96 of mouse CD59a (NP_031678.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O55186
Protein Weight	14kDa
Shipping	Ice bag