

Rabbit Anti-Mouse Cd6 Monoclonal Antibody For Fc - Q61003

Item Number: CM0011111

Introduction

ABflo 647-conjugated rabbit anti-mouse CD6 monoclonal antibody validated for flow cytometry (FC) in mouse samples. Targets T-cell differentiation antigen CD6, involved in T-cell activation and immune response. Suitable for research use in immunology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD6 mAb
Remarks/Alias	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Species	Mouse
GeneID (human)	923
GeneID	12511
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-155 of mouse CD6 (NP_033982.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61003
Protein Weight	64kDa/68kDa/72kDa
Shipping	Ice bag