

Rabbit Anti-Mouse Cd6 Monoclonal Antibody For Flow Cytometry - Q61003

Item Number: CM0009984

Introduction

-conjugated rabbit monoclonal antibody against mouse CD6 for flow cytometry. Recombinant immunogen (aa 20-155). Useful for investigating T-cell costimulation, immunological synapse maturation, and bacterial pattern recognition. - Q61003

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD6 mAb
Remarks/Alias	Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Myomorpha; Muroidea; Muridae; Murinae; Mus; Mus
Species	Mouse
GeneID (human)	923
GeneID	12511
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-155 of mouse CD6(NP_033982.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q61003
Protein Weight	64kDa/68kDa/72kDa
Shipping	Ice bag