

Rabbit Anti-Mouse Cd62L Monoclonal Antibody For Flow Cytometry - P18337

Item Number: CM0008867

Introduction

High-quality rabbit anti-mouse CD62L (L-selectin) monoclonal antibody, validated for flow cytometry. Suitable for lymphocyte homing, leukocyte adhesion, and immune cell trafficking studies. Recombinant immunogen; recognizes mouse CD62L; target 42 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD62L mAb
Remarks/Alias	Lnhr; CD62L; LAM-1; Ly-22; Lyam1; LECAM1; Ly-m22; Lyam-1; LECAM-1; L-selectin
Species	Mouse
GeneID (Human)	6402
GeneID	20343
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P18337
Protein Weight	42kDa
Shipping	Ice bag