

Rabbit Anti-Mouse Cd68 Mab For If/Icc And Fc - P31996

Item Number: CM0012365

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD68 (macrosialin) for IF/ICC and FC applications. Recognizes a lysosomal/endosomal membrane protein expressed in macrophages and dendritic cells. - P31996.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD68 mAb
Remarks/Alias	Lamp4; gp110; Scard1
Species	Mouse
GeneID (Human)	968
GeneID	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	34KDa
Shipping	Ice bag