

Rabbit Anti-Mouse Cd68 Monoclonal Antibody For Flow Cytometry - P31996

Item Number: CM0012623

Introduction

Rabbit anti-Mouse CD68 monoclonal antibody validated for intracellular flow cytometry in mouse samples. Targets macrophage marker CD68 (P31996). Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD68 mAb
Remarks/Alias	Lamp4; gp110; Scard1
Species	Mouse
GeneID (Human)	968
GeneID (Mouse)	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	35kDa
Shipping	Ice bag