

Rabbit Anti-Mouse Cd68 Monoclonal Antibody For If/Icc, Fc - P31996

Item Number: CM0011734

Introduction

High-specificity -conjugated rabbit monoclonal antibody targeting mouse macrosialin (CD68) for immunofluorescence and flow cytometry. This dye-labeled antibody provides reliable detection of CD68 expression in macrophages and dendritic cells, supporting immunological studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD68 mAb
Remarks/Alias	Lamp4; gp110; Scard1
Species	Mouse
Gene ID (Human)	968
Gene ID	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	33KDa/34KDa
Shipping	Ice bag