

Rabbit Anti-Mouse Cd68 Monoclonal Antibody For Intracellular Flow Cytometry - P31996

Item Number: CM0012121

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD68 (macrosialin, Lamp4, gp110). Validated for intracellular flow cytometry. Reacts with mouse. Ideal for studying macrophage and dendritic cell phagocytic activities.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD68 mAb
Remarks/Alias	Lamp4; gp110; Scard1
Species	Mouse
GeneID (Human)	968
GeneID	12514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-291 of mouse CD68 (NP_001277987.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Mouse
SWISS	P31996
Protein Weight	35kDa
Shipping	Ice bag