

Rabbit Anti-Mouse Cd69 Mab For Fc - P37217

Item Number: CM0004011

Introduction

-conjugated rabbit monoclonal antibody against mouse CD69, suitable for flow cytometry (FC). Validated for mouse reactivity; ideal for studying early T-cell activation and immune cell homeostasis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
GeneID (Human)	969
GeneID	12515
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to aminoacids 64-199 of mouse CD69(NP_001028294.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22kDa
Shipping	Ice bag