

Rabbit Anti-Mouse Cd69 Monoclonal Antibody For Fc - P37217

Item Number: CM0010542

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD69 (AIM/VEA) for flow cytometry applications. Reacts with mouse samples. Ideal for T-cell activation and immune cell homeostasis studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
GeneID (Human)	969
GeneID	12515
Immunogen	Recombinant Protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22kDa
Shipping	Ice bag