

Rabbit Anti-Mouse Cd69 Monoclonal Antibody For Fc - P37217

Item Number: CM0004753

Introduction

rabbit monoclonal antibody targeting mouse CD69 (AIM/VEA), validated for flow cytometry. Detects early activation antigen on T-cells, NK cells, and other immune cell subsets, suitable for immuno-oncology and inflammation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD69 mAb
Remarks/Alias	AIM; VEA; 5830438K24Rik
Species	Mouse
GeneID (Human)	969
GeneID	12515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 64-199 of mouse CD69 (NP_001028294.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P37217
Protein Weight	22 kDa
Shipping	Ice bag