

Rabbit Anti-Mouse Cd71 Monoclonal Antibody For Flow Cytometry - Q62351

Item Number: CM0012774

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD71 (transferrin receptor 1). Validated for flow cytometry. Ideal for research on iron metabolism, cell proliferation, and transferrin receptor trafficking.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD71 mAb
Remarks/Alias	TR, TFR, p90, CD71, TFR1, Trfr, Mtvr1, 2610028K12Rik, E430033M20Rik
Species	Mouse
GeneID (Human)	7037
GeneID	22042
Immunogen	Recombinant protein corresponding to a sequence within amino acids 89-763 of mouse CD71 (NP_035768.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q62351
Protein Weight	86kDa
Shipping	Ice bag