


Rabbit Anti-Mouse Cd72 Mab For Fc - P21855

Item Number: CM0012679



Introduction

-conjugated rabbit monoclonal antibody against mouse CD72 (Lyb-2/Ly-32), validated for flow cytometry. Ideal for investigating B-cell receptor co-receptor functions, regulatory T/B-cell homeostasis, and B-cell differentiation in mouse research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD72 mAb
Remarks/Alias	CD72c; Ly-19; Ly-32; Lyb-2; Ly-m19
Species	Mouse
GeneID (Human)	971
GeneID	12517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 226-354 of mouse CD72 (NP_031680.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P21855
Protein Weight	40kDa
Shipping	Ice bag