


Rabbit Anti-Mouse Cd72 Mab - P21855

Item Number: CM0004116



Introduction

Rabbit anti-Mouse CD72 mAb for flow cytometry. Recognizes mouse CD72 (Lyb-2, Ly-32), a B-cell differentiation antigen and BCR co-receptor. Monoclonal antibody sourced from rabbit, ideal for B-cell immunophenotyping in research settings.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD72 mAb
Remarks/Alias	CD72c; Ly-19; Ly-32; Lyb-2; Ly-m19
Species	Mouse
GeneID (Human)	971
GeneID	12517
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P21855
Protein Weight	40kDa
Shipping	Ice bag