


Rabbit Anti-Mouse Cd72 Monoclonal Antibody For Fc - P21855

Item Number: CM0011174



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD72 for flow cytometry. Validated for mouse cross-reactivity, molecular weight 40kDa, suitable for B-cell co-receptor and B-cell regulation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD72 mAb
Remarks/Alias	CD72c; Ly-19; Ly-32; Lyb-2; Ly-m19
Species	Mouse
GeneID (human)	971
GeneID	12517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 226-354 of mouse CD72 (NP_031680.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P21855
Protein Weight	40kDa
Shipping	Ice bag