

Rabbit Anti-Mouse Cd79B Monoclonal Antibody For Flow Cytometry - P15530

Item Number: CM0009453



Introduction

-conjugated rabbit monoclonal antibody against mouse CD79b (B29/Ig-beta) validated for flow cytometry. Ideal for B-cell receptor signaling studies and immune phenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD79b mAb
Remarks/Alias	B29; Igb; Igbeta; Ig-beta
Species	Mouse
GeneID (Human)	974
GeneID	15985
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-158 of mouse CD79b (NP_032365.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P15530
Protein Weight	25kDa
Shipping	Ice bag