

Rabbit Anti-Mouse Cd80 Monoclonal Antibody For Flow Cytometry - Q00609

Item Number: CM0009216



Introduction

Rabbit anti-Mouse CD80 Monoclonal Antibody validated for flow cytometry. Detects mouse CD80 (B7-1), a costimulatory molecule on activated B-cells and monocytes. Suitable for T-cell activation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD80 mAb
Remarks/Alias	B71; Ly53; TSA1; Cd28l; Ly-53; MIC17
Species	Mouse
GeneID (human)	941
GeneID	12519
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-245 of mouse CD80 (NP_033985.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00609
Protein Weight	34kDa/24kDa
Shipping	Ice bag