

Rabbit Anti-Mouse Cd80 Monoclonal Antibody For Flow Cytometry - Q00609

Item Number: CM0009896



Introduction

Rabbit anti-Mouse CD80 Monoclonal Antibody (FC) targets mouse CD80/B7-1, a key costimulatory molecule in T-cell activation. Recombinant immunogen ensures specificity. Ideal for immune cell analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD80 mAb
Remarks/Alias	B71; Ly53; TSA1; Cd28l; Ly-53; MIC17
Species	Mouse
GeneID (Human)	941
GeneID	12519
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 38-245 of mouse CD80(NP_033985.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q00609
Protein Weight	34kDa
Shipping	Ice bag