

Rabbit Anti-Mouse Cd82 Monoclonal Antibody For Flow Cytometry - P40237

Item Number: CM0009250

Introduction

-conjugated rabbit monoclonal antibody against mouse CD82 (C33/IA4/Kai1/Tspan27), validated for flow cytometry. Target is a tetraspanin involved in signal transduction, adhesion, and migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD82 mAb
Remarks/Alias	C33; IA4; Kai1; Tspan27
Species	Mouse
GenelD (Human)	3732
GenelD	12521
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-266 of mouse CD82 (NP_001258390.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P40237
Protein Weight	30kDa
Shipping	Ice bag