

## Rabbit Anti-Mouse Cd83 Monoclonal Antibody For Fc - O88324

Item Number: CM0005959

### Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD83, validated for flow cytometry. Ideal for dendritic cell profiling and immune response studies.

[Learn More](#)

| Parameter        | Value                                                                                                             |
|------------------|-------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD83 mAb                                                                                        |
| Remarks/Alias    | mCD83                                                                                                             |
| Species          | Mouse                                                                                                             |
| Gene ID (Human)  | 9308                                                                                                              |
| Gene ID          | 12522                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 22-133 of mouse CD83 (NP_033986.1). |
| Source           | Rabbit                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                             |
| Application      | FC                                                                                                                |
| Cross Reactivity | Mouse                                                                                                             |
| SWISS            | O88324                                                                                                            |
| Protein Weight   | 21KDa                                                                                                             |
| Shipping         | Ice bag                                                                                                           |