

Rabbit Anti-Mouse Cd83 Monoclonal Antibody For Flow Cytometry - O88324

Item Number: CM0012080

Introduction

-conjugated rabbit monoclonal antibody against mouse CD83, validated for flow cytometry. This antibody targets the extracellular domain and is suitable for immune cell profiling and CD83 expression studies in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD83 mAb
Remarks/Alias	mCD83
Species	Mouse
Gene ID (Human)	9308
Gene ID	12522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-133 of mouse CD83 (NP_033986.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88324
Protein Weight	21KDa
Shipping	Ice bag