

Rabbit Anti-Mouse Cd83 Monoclonal Antibody For Flow Cytometry - O88324

Item Number: CM0012300

Introduction

-conjugated rabbit monoclonal antibody against mouse CD83, validated for flow cytometry. Detects CD83, a transmembrane glycoprotein involved in immune regulation, expressed on dendritic cells and lymphocytes. For research use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD83 mAb
Remarks/Alias	mCD83
Species	Mouse
GeneID (Human)	9308
GeneID	12522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-133 of mouse CD83 (NP_033986.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88324
Protein Weight	21KDa
Shipping	Ice bag