

Rabbit Anti-Mouse Cd83 Monoclonal Antibody For Flow Cytometry - O88324

Item Number: CM0011987

Introduction

Rabbit anti-Mouse CD83 monoclonal antibody validated for flow cytometry, targeting mouse CD83 antigen (O88324). Specific binding to mouse immune cells for research use.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD83 mAb
Remarks/Alias	mCD83
Species	Mouse
Gene ID (Human)	9308
Gene ID	12522
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-133 of mouse CD83 (NP_033986.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	O88324
Protein Weight	21KDa
Shipping	Ice bag