


Rabbit Anti-Mouse Cd87 Monoclonal Antibody For Flow Cytometry - P35456

Item Number: CM0012488



Introduction

-conjugated rabbit monoclonal antibody against mouse CD87 (uPAR) for flow cytometry. Validated for mouse samples. Targets urokinase receptor involved in cell migration and tissue remodeling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD87 mAb
Remarks/Alias	Cd87; uPAR; u-PAR
Species	Mouse
GeneID (human)	5329
GeneID	18793
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-297 of mouse CD87 (NP_035243.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35456
Protein Weight	24kDa/35kDa
Shipping	Ice bag