

Rabbit Anti-Mouse Cd88/C5Ar Mab - P30993

Item Number: CM0007761



Introduction

-conjugated rabbit monoclonal antibody against mouse CD88/C5aR (C5aR1) validated for flow cytometry. Targets the N-terminal domain of this G protein-coupled receptor involved in inflammation and immune chemotaxis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD88/C5aR mAb
Remarks/Alias	C5aR; C5r1; Cd88; D7Msu1
Species	Mouse
Gene ID (Human)	728
Gene ID	12273
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-37 of mouse CD88/C5aR (NP_031603.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P30993
Protein Weight	39kDa
Shipping	Ice bag