

Rabbit Anti-Mouse Cd88/C5Ar Monoclonal Antibody For Flow Cytometry - P30993

Item Number: CM0011760



Introduction

-conjugated rabbit anti-mouse CD88/C5aR monoclonal antibody for flow cytometry. Detects 39kDa C5a anaphylatoxin chemotactic receptor 1 (P30993) in mouse samples. Ideal for inflammation and chemotaxis research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD88/C5aR mAb
Remarks/Alias	C5aR; C5r1; Cd88; D7Msu1
Species	Mouse
GeneID (Human)	728
GeneID	12273
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-37 of mouse CD88/C5aR (NP_031603.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P30993
Protein Weight	39kDa
Shipping	Ice bag