


Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Fc - P01731

Item Number: CM0011936



Introduction

High-quality Rabbit anti-Mouse CD8a monoclonal antibody for flow cytometry (FC). Reacts specifically with mouse CD8a (Ly-2, Ly-35; SWISS: P01731). Research-grade for immune profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
GeneID (Mouse)	12525
GeneID (Human)	925
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag