

## Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Flow Cytometry - P01731

Item Number: CM0012639



### Introduction

-conjugated rabbit monoclonal antibody against mouse CD8a (Lyt-2), validated for flow cytometry (FC). Ideal for immuno-phenotyping of mouse T lymphocytes in research and diagnostic applications.

[Learn More](#)

| Parameter        | Value                                                                                                                                    |
|------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Mouse CD8a mAb                                                                                                               |
| Remarks/Alias    | Ly-2; Ly-B; Ly-35; Lyt-2                                                                                                                 |
| Species          | Mouse                                                                                                                                    |
| Gene ID (Human)  | 925                                                                                                                                      |
| Gene ID          | 12525                                                                                                                                    |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1). |
| Source           | Rabbit                                                                                                                                   |
| Category         | Monoclonal Antibodies                                                                                                                    |
| Application      | FC                                                                                                                                       |
| Cross Reactivity | Mouse                                                                                                                                    |
| SWISS            | P01731                                                                                                                                   |
| Protein Weight   | 27KDa                                                                                                                                    |
| Shipping         | Ice bag                                                                                                                                  |