

Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Flow Cytometry - P01731

Item Number: CM0010327



Introduction

Rabbit anti-Mouse CD8a Monoclonal Antibody for flow cytometry. Recognizes T-cell surface glycoprotein CD8 alpha chain (P01731), suitable for mouse immune response studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
Gene ID (Human)	925
Gene ID	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27KDa
Shipping	Ice bag