


Rabbit Anti-Mouse Cd8A Monoclonal Antibody For Flow Cytometry - P01731

Item Number: CM0010013



Introduction

High-quality -conjugated rabbit monoclonal antibody against mouse CD8a, validated for flow cytometry. Ideal for immunology research on T-cell subsets and cytotoxic responses. Reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD8a mAb
Remarks/Alias	Ly-2; Ly-B; Ly-35; Lyt-2
Species	Mouse
Gene ID (Human)	925
Gene ID (Mouse)	12525
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-196 of mouse CD8a (NP_001074579.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01731
Protein Weight	27kDa
Shipping	Ice bag