

Rabbit Anti-Mouse Cd90.2 Monoclonal Antibody For Fc - P01831

Item Number: CM0009112

Introduction

-labeled rabbit monoclonal antibody against mouse CD90.2 (Thy-1), validated for flow cytometry (FC). Ideal for immunophenotyping and studying cell-cell interactions in mouse immune and neuronal systems.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD90.2 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
GeneID (Human)	7070
GeneID	21838
Immunogen	Recombinant protein of mouse CD90.2
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag