

Rabbit Anti-Mouse Cd90.2 Monoclonal Antibody For Flow Cytometry - P01831

Item Number: CM0010359

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD90.2 (Thy-1.2). Ideal for sensitive and specific flow cytometry detection of mouse Thy-1 glycoprotein expression in immunological studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD90.2 mAb
Remarks/Alias	T25; CD90; Thy-1; Thy1.1; Thy1.2; Thy-1.2
Species	Mouse
GeneID (Human)	7070
GeneID	21838
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P01831
Protein Weight	18kDa
Shipping	Ice bag