


Rabbit Anti-Mouse Cd95/Fas Monoclonal Antibody For Flow Cytometry - P25446

Item Number: CM0006077



Introduction

-conjugated rabbit monoclonal antibody against mouse CD95/FAS (P25446) for flow cytometry. Detects mouse FAS receptor, involved in apoptosis and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD95/FAS mAb
Remarks/Alias	lpr; APO1; APT1; CD95; TNFR6; Tnfrsf6
Species	Mouse
GeneID (Human)	355
GeneID	14102
Immunogen	Recombinant protein of mouse CD95/FAS.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P25446
Protein Weight	37KDa
Shipping	Ice bag