

Rabbit Anti-Mouse Cd97 Monoclonal Antibody For Flow Cytometry - Q9Z0M6

Item Number: CM0011777

Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD97 (Adgre5), validated for flow cytometry. Suitable for studying leukocyte adhesion, migration, and immune cell profiling in mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD97 mAb
Remarks/Alias	Cd97; TM7LN1
Species	Mouse
GeneID (Human)	976
GeneID (Mouse)	26364
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-265 of mouse CD97 (NP_036055.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9Z0M6
Protein Weight	90KDa
Shipping	Ice bag