

Rabbit Anti-Mouse Cd98 Monoclonal Antibody For Flow Cytometry - P10852

Item Number: CM0008158



Introduction

-conjugated rabbit monoclonal antibody targeting mouse CD98 (SLC3A2/4F2), validated for flow cytometry. Binds the 58 kDa amino acid transporter heavy chain, cross-reacts with mouse. Suitable for immunophenotyping and transporter expression studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CD98 mAb
Remarks/Alias	4F2; Cd98; Ly10; Mdu1; 4F2HC; Ly-10; NACAE; Ly-m10; Mgp-2hc
Species	Mouse
GeneID (Human)	6520
GeneID (Mouse)	17254
Immunogen	Recombinant protein corresponding to amino acids 105-526 of mouse CD98 (NP_032603.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Mouse
SWISS	P10852
Protein Weight	58kDa
Shipping	Ice bag