

Rabbit Anti-Mouse Cx3Cr1 Monoclonal Antibody For Fc - Q9Z0D9

Item Number: CM0011438

Introduction

-conjugated rabbit monoclonal antibody against mouse CX3CR1, validated for flow cytometry. Suitable for investigating fractalkine receptor roles in immune cell migration, microglial function, and gut immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse CX3CR1 mAb
Remarks/Alias	mCX3CR1
Species	Mouse
Gene ID	13051
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse CX3CR1 (NP_034117.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9Z0D9
Protein Weight	40 kDa
Shipping	Ice bag