

Rabbit Anti-Mouse Cxcr2/Il-8Rb Monoclonal Antibody For Fc - P35343

Item Number: CM0009413

Introduction

-conjugated rabbit anti-mouse CXCR2/IL-8RB monoclonal antibody for flow cytometry. Targets the chemokine receptor CXCR2, crucial for neutrophil trafficking and inflammation studies. Ice bag shipping.

[Learn More](#)

Specification	Details
Product Name	Rabbit anti-Mouse CXCR2/IL-8RB mAb
Remarks/Alias	CD128; IL8RA; IL8rb; CDw128; Cmkar2; Gpcr16; IL-8Rh; IL-8rb; mL-8RH
Species	Mouse
Gene ID (Human)	3579
Gene ID	12765
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-44 of mouse CXCR2/IL-8RB (NP_034039.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P35343
Protein Weight	40kDa
Shipping	Ice bag