


Rabbit Anti-Mouse E-Cadherin/Cd324 Monoclonal Antibody For Fc - P09803

Item Number: CM0006034



Introduction

-conjugated rabbit monoclonal antibody targeting mouse E-Cadherin/CD324, validated for flow cytometry. Specifically binds mouse E-Cadherin with no cross-reactivity reported. Ideal for epithelial cell adhesion and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse E-Cadherin/CD324 mAb
Remarks/Alias	Um; UVO; Ecad; ARC-1; E-cad; L-CAM
Species	Mouse
GeneID (Human)	999
GeneID	12550
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 157-369 of mouse E-Cadherin/CD324 (NP_033994.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P09803
Protein Weight	98kDa
Shipping	Ice bag