

Rabbit Anti-Mouse FcγRI/CD64 Mab - P26151

Item Number: CM0012575

Introduction

High-affinity rabbit monoclonal antibody targeting mouse FcγRI/CD64, validated for flow cytometry. Detects mouse CD64, a key receptor in IgG-mediated innate and adaptive immunity. Suitable for mouse macrophage studies, this antibody enables precise immune profiling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Mouse FcγRI/CD64 mAb
Remarks/Alias	CD64; IGGHAFC; FcγgammaRI
Species	Mouse
GeneID (Human)	2209
GeneID (Mouse)	14129
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 25-297 of mouse FcγRI/CD64 (NP_034316.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	P26151
Protein Weight	44kDa
Shipping	Ice bag